Vital Base Sit

Vital Base Royal

Vital Base Venti

Deutsch	Bedienungsanleitung
English	Operation Instructions
Italiano	Istruzioni d'uso
Français	Mode d'emploi
Norsk	Bruerveidningen
Dansk	Brugsanvisning

Deutsch

INHALT

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English

CONTENTS

1. Purpose of product
 2. Range of use
 3. Technical specifications
 4. Potential risk
 5. Mounting description
 6. Maintenance
 7. Routines for inspection
 8. Washing/Reuse
-

1. Purpose of product

VB Sit, VB Royal and VB Venti are seat cushions developed to fit all types of wheelchairs. They are built on the concept of simplicity in use in combination with good pressure ulcer preventive properties. The pressure ulcer preventive properties cover both pressure distribution and micro climate.

A combination of a friction textile and velcro under the cushion makes it compatible with all existing wheelchairs on the market.

2. Range of use

VB Sit is a pressure ulcer preventive cushion built from different layers of polyurethane foam for long time pressure distribution. It has side supports for increased stability and

contact area. An integrated abduction block offers a good foundation for erection of the upper trunk. Integrated support area in front of the seat bones.

The cushion is well suited for pressure ulcer prevention related to long time seating, and where one is seeking stability for the user.

VB Royal is a pressure ulcer preventive cushion built from different layers of polyurethane foam for long time pressure distribution. An integrated PU-Gel is used for the purpose of shock absorption. An integrated abduction block offers a good foundation for erection of the upper trunk. Integrated support area in front of the seat bones.

Outer cover with temperature stabilizing properties supplied as standard.

The cushion is well suited for pressure ulcer prevention related to long time seating, and where one is seeking pelvic stability and shock absorbing properties for the user.

VB Venti is a pressure ulcer preventive cushion with extremely high ventilating properties. VB Venti is designed for easy usage. It is built from three layers of mesh fabric, where both sides provides the same good pressure distributive properties. This reduces possibilities for incorrect usage.

3. Technical specifications

VB Sit

CORE

4 layers polyurethane foam. The cushion is made from polyurethane with different weight carrying properties. Two of the layers are mounted next to the other horizontally in the middle section. At the rear is used a softer polyurethane foam than in the front. This is to assure stabilization in front of the seat bones.

It has side supports for increased stability and contact area and an integrated abduction block for assistance for erection of the upper trunk.

INNER COVER

The core has a white protective cover, TechmaFlex, protecting the core from pollution. TechmaFlex is a semi permeable "breathing" textile made from polyurethane. It lets vapour through, while water and other liquids are prevented from entering the core. TechmaFlex is very elastic.

OUTER COVER

100% knitted polyester with 3D-structure. The cover is very durable and stretchable with a high degree of air penetration.

VB Royal

CORE

6 layers polyurethane foam + 1 polyurethane gel. The cushion is made from polyurethane with different weight carrying properties. 4 of the layers are placed horizontally next to

each other in different levels. At the rear is used a softer polyurethane foam than in the front. This is to assure stabilization in front of the seat bones. The PU-Gel assures shock absorption. I.

INNER COVER

The core has a white protective cover, TechmaFlex, protecting the core from pollution. TechmaFlex is a semi permeable "breathing" textile made from polyurethane. It lets vapour through, while water and other liquids are prevented from entering the core. TechmaFlex is very elastic.

OUTER COVER

Made from 3 layers. On top is a 100% knitted polyester. In the middle is a 3D Mesh for ventilation. At the bottom is Outlast, a PCM with temperature stabilizing properties. The cover is stretchable and durable, holds a high degree of ventilation and temperature stabilizing properties.

VB Venti

CORE

3 layers of 3D Mesh. 3D Mesh is a textile with very open structure. The textile has good elasticity and high ventilating properties. The core is built from 3 layers of 3D Mesh. These layers are encapsulated by a thinner 3D textile. This assures a good pressure distribution throughout the cushion. The open structure makes the cushion easy to maintain. The core can be washed in the shower, shaken and set standing to dry. The core dries fast.

OUTER COVER

100% knitted polyester with 3D-structure. The cover is very durable and stretchable with a high degree of air penetration.

The bottom textile of the outer covers are made from a friction fabric with mounted velcro, securing the cushion to stay in position in any wheelchair model and on any surface.

4. Potential risk

The cushions are "pre-fit", and needs no more fitting as long as mounting- and washing instructions are followed. If the user has open ulcers, we recommend the user not to sit on the cushion.

5. Mounting description

The cushion is taken out of the packing and placed in the seat of the wheelchair. Place the cushion according to instructions on the label attached to the bottom fabric of the outer cover. Check that the velcro and friction fabric is facing down towards the seat plate. The front of the cushion should be even with the front of the seat plate.

6. Maintenance

The only need for maintenance of the cushions is cleaning:

1. Open the zipper at the rear and take the outer cover off.

2. Open the zipper at the bottom of the middle white TechmaFlex cover (VB Sit & VB Royal).
3. Take out the core, and place it edgewise for airing while the covers are washed.
4. Follow washing instructions described on label attached to the bottom of the middle and outer cover.
5. The core is placed in the covers with the label facing down. For correct position of front and rear, see label.
6. Make sure that the label and zipper of the middle cover (if used) is facing down.
7. Mount the outer cover. Make sure that the zipper faces rearwards and that the label faces down.

7. Routines for inspection

Before seating the user, check that the cushion is placed correct according to instruction on label under the cushion. On a regular basis, check that the core and covers are placed correctly (follow description under chapters 5 & 6).

The need for washing core & covers will vary from user to user depending on degree of incontinence etc. A rule of thumb can be to take the core out of the covers once a week and let it stand edgewise to aerate and at the same time wash the covers.

It is recommended to have an extra cover for the cushion so that the user doesn't need to wait for the covers to dry.

The core and the covers are to be washed seperately. Use mild detergent. DO NOT USE

CHLORINE BLEACH. After washing, shake the core to remove excess water and let it stand to dry. Do not use harsh cleaners such as scouring pads & heavy-duty grease removers,

or solvents, such as acetone. To remove spots and stains, use standard household cleansers and a soft cloth. To loosen heavy, dried-on soil, you may first need to saturate the spot.

8. Washing / Reuse

	VB Sit	VB Royal	VB Venti
CORE			
Wash	Hand wash 40 °C	Hand wash 40 °C	Hand wash 40 °C
Disinfection	Virkon S*	Virkon S*	Virkon S*
	Auto clave 105 °C	Auto clave 105 °C	Auto clave 105 °C
Drying	Squeeze	Squeeze	Squeeze
	Air dry standing edgewise	Air dry standing edgewise	Air dry standing edgewise
INNER COVER			Not standard
Wash	Machine wash 95 °C	Machine wash 95 °C	
Drying	Tumble dry max 85 °C	Tumble dry max 85 °C	
OUTER COVER			
Wash	Machine wash 60 °C	Machine wash 60 °C	Machine wash 60 °C
Drying	Tumble dry max 85 °C	Tumble dry max 85 °C	Tumble dry max 85 °C

* Read instructions on package.

Alu Rehab considers seat cushions to be an intimate and very personal technical aid. If reusing, we advice to at least exchange the outer and also the inner covers.

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