

Product Info

Date: 20/6 - 2006	Product info.: 7-2006	Written by: Fredrik Delås
Item: New arm support pads	Delivery: From this date	

From this date we will start to deliver new arm supports.



The arm support has gone through a "face lift" and is now made of moulded polyurethane. It is easy to clean, and has no plastic edges or cracks where dirt has a tendency to find its way. The first arm support pads of this model are mounted in production now.

Item number	Item
20835	Arm support pad complete with screws
21682	Arm support left complete
21681	Arm support right complete