

Product Info

| | | |
|-------------------------|-----------------------|---------------------------|
| Date: 31/5 - 2007 | Product info.: 2-2007 | Written by: Fredrik Delås |
| Item: Lille Viking 2000 | | |

There is now signed a combination agreement between U-B-Let A/S & Alu Rehab A/S.

U-B-Let A/S assistive engine Lille Viking 2000 can now be used on Netti I, Netti II, Netti III, Netti Svipp, Netti comfort 4U|CE, Netti comfort 4U|CES & Netti comfort 4U|CED. There must be used two different models of Lille Viking 2000 on the different Netti chairs:

Netti I, Netti II, Netti III & Netti Svipp: Lille Viking 2000 LOW modell.
Netti 4U comfort|CE/CES/CED: Lille Viking 2000 FLAT modell.

Mounting descriptions are sent separately from the producer U-B-Let A/S.