

PRODUCT INFORMATION

Date: 28/1 - 2008	Product info.: 3-2008	Written by: Fredrik Delås
Item: Friction fabric on cushions	Delivery: From this date	

All Vital Base and Alu Rehab seat cushions will now be supplied with a friction fabric

The friction fabric will be used as the bottom textile of all seat cushions. It has a surface which provides excellent stabilization of the cushion against the seat. This replaces the need for velcro, and makes the cushions much more versatile for usage in other wheelchairs. The change starts from this date, but there will in a period of time be delivered cushions with velcro. This depends on numbers and sizes of cushions with velcro left on stock.

The friction fabric fulfills the same documentation and demands as the existing outer covers.

Due to that our wheelchairs and cushions are being re-used in most markets, we will still produce seat plates and sling seats with velcro until further notice.

This change does NOT imply change of item numbers.