

# BORD - NETTI RULLESTOL

## Hvilket bord passer til hvilken rullestolmodell?

		PASSER		PASSER IKKE		RULLESTOLENE ER MED STANDARD ARMLEN										
PASSER HVILKEN RULLESTOL OG SETEBREDE	BORD															
	Netti Mini	Universal lite/mini - art.nr 15102	Universal medium - art.nr 15100	Universal stort - art.nr 15101	Universal lås nedfellbart - art.nr - 15139 lite/mini	Universal nedfellbart - art.nr 15140 medium	Universal nedfellbart - art.nr 15141 stort	Universal låsbart lite/mini - art.nr 15123	Universal låsbart medium - art.nr 15124	Universal låsbart stort - art.nr 15125	Plexi bord 55 cm - art.nr 21608	Plexi bord 60 cm - art.nr 21609	Universal nedfellbart CED/ CEDs medium - art 28826	Universal nedfellbart stort CED/CEDs - art 28832	Universal låsbart medium CED/CEDs - artnr 28838	Universal låsbart stort CED/CEDs - artnr 28836
	Netti Mini	Netti Mini 25, 30 cm	Netti Mini 35 cm	Netti Mini	Netti Mini 25, 30 cm	Netti Mini 35 cm	Netti Mini	Netti Mini 25, 30 cm	Netti Mini 35 cm	Netti Mini						
	Netti I/II/III	Netti I/II/III 35, 38, 40 cm	Netti I/II/III 43, 45, 50 cm	Netti I/II/III	Netti I/II/III	Netti I/II/III 35, 38, 40 cm	Netti I/II/III 43, 45, 50 cm	Netti I/II/III	Netti I/II/III 35, 38, 40 cm	Netti I/II/III 43, 45, 50 cm	Netti I/II/III 55 cm with wide arm support	Netti I/II/III 60 cm with wide arm support				
	Netti 4U CE/CES/FA	Netti 4U CE/CES/FA 35, 40 cm	Netti 4U CE/CES/FA 45, 50 cm	Netti 4U CE/CES/FA	Netti 4U CE/CES/FA 35, 40 cm	Netti 4U CE/CES/FA 45, 50 cm	Netti 4U CE/CES/FA	Netti 4U CE/CES/FA 35, 40 cm	Netti 4U CE/CES/FA 45, 50 cm	Netti 4U CE/CES/FA 45, 50 cm						
	Netti 4U CED/CEDS	Netti 4U CED/CEDS	Netti 4U CED/CEDS	Netti 4U CED/CEDS	Netti 4U CED/CEDS	Netti 4U CED/CEDS	Netti 4U CED/CEDS	Netti 4U CED/CEDS	Netti 4U CED/CEDS	Netti 4U CED/CEDS			Netti 4U CED/CEDS 35, 40 cm	Netti 4U CED/CEDS 43, 45 50 cm	Netti 4U CED/CEDS 35, 40 cm	Netti 4U CED/CEDS 43, 45 50 cm
Netti Svipp	Netti Svipp 35, 38 40 cm	Netti Svipp 43, 45 50 cm	Netti Svipp	Netti Svipp 35, 38 40 cm	Netti Svipp 43, 45 50 cm	Netti Svipp	Netti Svipp 35, 38 40 cm	Netti Svipp 43, 45 50 cm	Netti Svipp							



Universal bord kommer i 3 str.

- Navn: Universal bord mini  
for setebredde: 25-30 cm  
Art nr: 15102  
Navn: Universal bord medium  
for setebredde: 35 - 43 cm  
Art. nr: 15100  
Navn: Universal bord stort  
for setebredde: 43-50 cm  
Art. nr: 15101



Universal bord med lås

- Navn: Universal bord med lås, mini  
for setebredde: 25-35 cm  
Art. nr: 15123  
Navn: Universal bord med lås,  
medium, for sb: 35-40 cm  
Art. nr: 15124  
Navn: Universal bord med lås, stort  
for setebredde 43 -50 cm  
Art. nr: 15125



Universalbord nedfellbart og med lås

- Navn: Universalbord nedfellbart, mini  
for setebreddene: 25-35 cm  
Art.nr: 15139  
Navn: Universalbord nedfellbart  
Medium for sb: 35-40 cm  
Art nr: 15140  
Navn: Universalbord nedfellbart, stort  
for setebreddene: 43-50 cm  
Art. nr: 15141